



FALL DETECTION AND FALL PREVENTION SUPPORT

Unique systems to detect users falling, wandering, suffering panic attacks or risking hypothermia all in a discreet way

- Detects a user falling and lying unconscious
- “Knock-over alarm” to summon assistance or call for assistance when pressed
- Monitors bed and chair occupancy
- Sensors set for times of immediate or prolonged vacation
- Releases carers from constant supervision
- Suitable for carers at home or professional care organisations



Stress free caring...



Made in the UK



ALARM DETECTOR



MOTION DETECTOR



CHAIR OCCUPANCY MONITOR BED OCCUPANCY MONITOR



ADVANCED RADIO PAGER

Sense-it monitoring systems offer a range of support for carers depending on individual need. They are easy to install and simple and secure to operate. Alarms are transmitted through a failsafe radio link to an Advanced Radio Pager or to an existing call system.

Alarm Detector

A unique Personal Alarm Badge which can detect falling, hyperthermia risk, panic or the user wandering too far from the Pager.

Motion Detector

A discreet detection unit with an infra red motion detector remote from the bed. The alarm is raised as soon as the users foot descends to the floor. The same sensor can be positioned to detect them sitting up, leaving a room, or entering a dangerous area. It can also detect door usage, where a door may be left open and replacing door mats which can be a trip hazard.

Bed and Chair Occupancy Monitor

A discreet sensor that detects the weight of the person lying in bed, or sitting in a chair. An alarm is raised either immediately on vacation, or after a prolonged period of vacation.

A unique feature is that if reset whilst the user is out of bed or a chair, it will stop detecting alarms until the use returns later. The occupancy mat can be changed for various other styles of sensor to suit the user need, please ask for details.

Advanced Radio Pager

All systems come with an Advanced Radio Pager. The pager distinguishes between the different severities of each alarm and is powered by rechargeable batteries.

For professional use a 32 channel revision is available. The pager has a 450m range.

We are the foremost supplier of Plesiocare Systems.

For further information or to order the right system, please call us on 0845 217 9952, or visit our website at www.alert-it.co.uk



ITS Designs Ltd, Atherstone House, Merry Lees Industrial Estate, Desford, Leicester LE9 9FE.

Telephone: 0845 217 9951 Fax: 0845 217 9953

Email: sales@alert-it.co.uk Web: www.alert-it.co.uk

Manufacturers and Distributors of Care Alarms