



The Schampi Pillow

Give Your Bed an Additional Sitting Function



In places such as homes for the elderly it is quite a common occurrence that care recipients want to sit on the edge of the bed, but are unable to. The Schampi pillow has been designed to solve this issue. The pillow is placed behind the back of the care recipient, once he or she has been helped in to a sitting position.

The pillow is prevented from sliding backward by securing the side rail of the bed behind the care recipient. It provides a secure and comfortable positioning thanks to the contoured back rest.

The Schampi Pillow

The user's self-esteem

Being able to sit on the edge of the bed is important for an individual's self-esteem. Becoming a "prisoner" in the supine position can be debilitating and depressing. It can provide increased well-being and feel stimulating to sit on the edge of the bed and eat or read the paper without assistance.

Moreover, many care recipients who suffer from anxiety/respiratory distress experience relief when they can sit on the edge of the bed for a moment. This could be applicable to care recipients with lung disease, such as asthma and COPD, or heart failure.

The Schampi pillow also has positive effects on pressure ulcer prevention, as it increases positioning options. Variation is made possible/encouraged. One-sided body positioning (supine) is avoided.

Provides relief for care staff

The Schampi pillow has been developed with the care recipient in mind, but also for the convenience of the care staff by releasing resources and improving ergonomics. The Schampi pillow provides support that can be helpful in feeding, dressing / undressing and other care duties.

Areas of use

All types of hospitals, homes for the elderly, service apartments, special needs housing, palliative care units and rehabilitation units. Other areas such as care homes for disabled and handicapped individuals, community-based housing, home care, and maternity wards.

FACTS

Cover

The pillow is available with a cover made of PU material. The cover is soft, comfortable and easy to clean. The cover is easy to remove for cleaning purposes (Velcro tape, see picture below).



Item no.

3354-PU Low back with armrest, PU cover
3356-PU Low back without armrest, PU cover

Size

40x70x85 cm
40x50x38 cm

Colour

blue/beige
blue/beige

Recycling

The Schampi pillow are fully recyclable and should be sorted as combustible waste. No dismantling is required.



The Schampi pillow must only be used on beds with side rails.



Design protection no: 002208306



Please read the instruction manual before using the product.



Järven Health Care

Bromsvägen 3 | SE-891 60 Örnsköldsvik | SWEDEN
Telephone: +46 (0)660-22 28 50 | Telefax: +46 (0)660 - 824 10
www.jarven.se | healthcare@jarven.se