

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

S3 has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 has a seat height of 47 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.  
**Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45-50 cm**

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function  
**Design your wheelchair by clicking in the circles in this guide.**  
**Blue, underlined text is clickable to show info** or image.  
**BB = seatwidth, HH= backrest height**

**Original design**

**Alternative options**

**Other options**

[info](#)

**1) Chassis colour**

- Black** - Anodic Black
- Gray** - Gray metallic
- Clear



**2) Seat width (BB)**

[art.no.](#)

- 33 cm
- 36 cm
- 39 cm
- 42 cm
- 45 cm
- 50 cm

- Clear

Art.no include: chassis, seat, calf strap, backrest lock



**3) Back rest**

- Standard**  
523BBHH VER2
- Clear

Height (HH)

- 20 cm
- 25 cm
- 30 cm
- 35 cm
- 40 cm
- 45 cm

- Clear

[Info](#)

- Tapered**  
525BBHH VER2
- Curved back**  
524BBHH VER2

[Info](#)

- Other backrest**

[List of options](#) art. \_\_\_\_\_



**4) Rear axle**

- Carbon standard**  
33023BB
- Clear

- Carbon +15 mm**  
33023BC\*

\* seat width BB +1

Use this rear axle when choosing arm support at Accessories and/or increased seat height

- Other rear axle or seat height**

[Increase the seat height](#)

[List of options](#) art. \_\_\_\_\_



**5) Castor and fork**

- 120 mm**  
[art.no.](#)
- Clear

- 90 mm**  
[art.no.](#)
- 150 mm**  
[art.no.](#)

**6) Footrest**

- Footrest Titanium**  
36001BB
- Footrest U-shape small heel bar**  
3600134 (width 33)
- Clear

- Titanium, Extended, reinforced**  
36111BB
- Foldable footplates**  
[art.no.](#)
- Heel strap for footplates**  
[art.no.](#)
- Plastic plate for footrest**  
36010BB

- Other footrest**

[List of options](#) art. \_\_\_\_\_



**7) Brake**

- Brake S3**  
4630000
- Clear

- One arm brake, right**  
[art.no.](#)

- Other brake**

[List of options](#) art. \_\_\_\_\_



**8) Rear wheel - without pushrim**

- 24" standard**  
[art.no.](#)
- Clear

- 22" Standard**  
[art.no.](#) [Rear wheel sizes](#)
- 25" Spinergy X** black spokes  
[art.no.](#)
- 26" Spinergy X** black spokes  
[art.no.](#)

- Other rear wheel**

[List of options](#) art. \_\_\_\_\_



**Tire- choose tire**

- Schwalbe Right Run** 22",24",25",26"
- Clear

- Rugged** 24"
- Schwalbe **Marathon Plus** 24",25",26"

- PU** 24" (massive)



**9) Pushrim**

- Titanium**  
190DD11 (x2 / wheelchair)

- Paragrip**  
19051DDBT x2
- Tetragrip**  
19052DDBTx2

- Clear **DD = Wheel diameter**

- Maxgrip** 22", 24"(contains latex)  
1903060 x2 (60=24", 80=22")
- Aluminum**  
1904404 x2 (04=24", 05=25", 06=26")

## Accessories



### Side guard

- [Side guard S3](#) [Video](#)  
[art.no.](#)
- Clear



### Arm support

- [Arm support S3](#) [Video](#)  
[art.no.](#)



### Side guard Carbon

- [Side guard Carbon](#)  
[art.no.](#)



### Push handles

- [Height adjustable S3](#) [Video](#)  
4252060 VER2
- [Height adjustable S3 +10](#)  
4252063 VER2  
Max. userweight 100kg
- Clear



### Push handle Foldable

- [Foldable S3](#) [Video](#)  
4245000



### Push bar

- [Push bar S3](#)  
42400BB VER2  
Max. userweight 100kg



### Anti-tip


- [Anti-tip](#) std seat height  
4353000 [Video](#)

[Other seat height](#)



### Cushion Wedge cushion

- [2,5 cm](#)
- [5 cm](#)
- [572BB40](#)
- [575BB40](#)
- Clear
- [Wedge cushion](#)  
57710BB



### Spoke protectors

- [24" logo](#)
- [22" logo](#)
- 4500000
- 4500100
- [24" no logo](#)
- [22" no logo](#)
- 4500011
- 4500111
- Clear



### Other

- [Tube protection front](#)  
9900251 (2 pcs)
- [Tube protection rear](#)  
\_9900251 (1 pcs)
- [Crutch holder S3](#) \_\_\_\_\_ pcs
- 9900301 (x1/crutch)
- [Backrest wedge](#) \_\_\_\_\_ st
- 4710012

## Set-up adjustments



### Back rest angle

Note measure (A) for desired back rest angle.  
14-28 mm. Standard=20 mm

A= \_\_\_\_\_ mm



### Rear axle adjustment

Note measure (A) for desired balance.  
Standard=60 mm, at raised rear axle= 51 mm


A= \_\_\_\_\_ mm



### Footrest height **Footbow - Foot plates**

Note measure (A) and (B) for desired footrest.  
[Bambino foot plate with Feewheel - info](#)

A= \_\_\_\_\_ mm      B= \_\_\_\_\_ mm



### Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle.



NOTE! See [manual](#) for more information.

Extra