

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g):

The weight changes based on the choices made below.

S3 large has an adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 large has a seat height of 49,5 cm in the front and 45,5 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options. **S3 large is 50 mm longer and the seat height is 25 mm higher than S3.**
Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45-50 cm

Save the guide to the "desktop" and open it with **Adobe Reader**, for best function
Design your wheelchair by clicking in the circles in this guide.
Blue, underlined text is clickable to show info or image.
BB = seatwidth, HH= backrest height



Original design

Alternative options

Other options

[info](#)

1) Chassis colour

- Black** - Anodic Black
- Clear

2) Seat width (BB)

[art.no](#)
 Clear

- 39 cm
- 42 cm
- 45 cm
- 50 cm

Art.no include:
 chassis, seat, calf
 strap, backrest lock

3) Back rest

- Standard**
523BBHH VER2
- Clear

Height (HH)

- 20 cm
- 25 cm
- 30 cm
- 35 cm
- 40 cm
- 45 cm
- Clear

Tapered

525BBHH VER2

Curved back

524BBHH VER2

Other backrest

[List of options](#) art. _____

4) Rear axle

- Carbon standard**
33023BB
- Clear

Carbon +15 mm

33023BC*

* seat width BB +1

Use this rear axle when choosing
 arm support at Accessories and/or increased seat height

Other rear axle or seat height

[Increase the seat height](#)
[List of options](#) art. _____

5) Castor and fork

- 120 mm**
[art.no.](#)
- Clear

90 mm

[art.no.](#)

150 mm

[art.no.](#)

6) Footrest

- Titanium, extended, reinforced**
36111BB
- Clear

Foldable footplates

[art.no.](#)

Heel strap for footplates

[art.no.](#)

Plastic plate for footrest

36010BB

Other footrest

[List of options](#) art. _____

7) Brake

- Brake S3**
4630000
- Clear

One arm brake, right

[art.no.](#)

Other brake

[List of options](#) art. _____

8) Rear wheel - without pushrim

- 24" standard**
[art.no.](#)
- Clear

24" Spinergy X black spokes [Rear wheel sizes](#)

[art.no.](#)

25" Spinergy X black spokes

[art.no.](#)

26" Spinergy X black spokes

[art.no.](#)

Other rear wheel

[List of options](#) art. _____

Tire - choose tire

- Schwalbe Right Run** 24",25",26"
- Clear

Rugged 24"

Schwalbe **Marathon Plus** 24",25",26"

PU 24" (massive)

9) Pushrim

- Titanium**
190DD11 (x2 / wheelchair)
- Clear

DD = Wheel diameter

Paragrip

19051DDBT x2

Tetragrip

19052DDBT x2

Aluminum

1904404 x2 (04=24", 05=25", 06=26")

Maxgrip 24" (contains latex)
1903060 x2

Accessories



Side guard

- [Side guard S3](#) [Video](#)
[art.no.](#)
- Clear



Arm support

- [Arm support S3](#) [Video](#)
[art.no.](#)



Side guard Carbon

- [Side guard Carbon](#)
[art.no.](#)



Push handles

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3 +10](#)
4252063 VER2
Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push bar

- [Push bar S3](#)
42400BB VER2
Max. userweight 100kg



Anti-tip

- [Anti-tip](#) std seat height
4353080 [Video](#)

[Other seat height](#)



Cushion Wedge cushion

- [2,5 cm](#) [Wedge cushion](#)
47245BB 57710BB
- [5 cm](#)
575BB40
- Clear



Spoke protectors

- [24" logo](#) [22" logo](#)
4500000 4500100
- [24" no logo](#) [22" no logo](#)
4500011 4500111
- Clear



Other

- [Tube protection front](#)
9900251 (2 pcs)
- [Tube protection rear](#)
9900251 (1 pcs)
- [Crutch holder S3](#) _____ pcs
- 9900301 (x1/crutch)
- [Backrest wedge](#) _____ st
4710012

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=60 mm, at raised rear axle= 51 mm


A= _____ mm



Footrest height [Footrest - Foot plates](#)

Note measure (A) and (B) for desired footrest.
[Bambino foot plate with Feewheel - info](#)

A= _____ mm B= _____ mm



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle. Art.no.9000152

NOTE! See [manual](#) for more information.



Extra