

GENERAL CATALOGUE DECUBITUS THERAPY

Your retailer

1.	ABOUT US	5 – 8
2.	HEMECARE	9 – 16
3.	SYSTEMATIC CLINIC CARE	17 – 38
4.	DECUBITUS THERAPY SEATING CUSHIONS	39 – 56
5.	POSITIONING AIDS	57 – 64

CONTENTS



ABOUT US



Company information



Funke Medical AG - since 1988

The Funke Medical AG is a medium-sized business with more than 100 employees and a total production area of 17.800m² and is based exclusively in Germany, in a small town called Raesfeld in North Rhine-Westphalia to be exact!

We are DIN EN ISO 13485 certified, which authorises us to manufacture medical devices and to sell them in accordance with relevant legislation.

Our quality management system is constantly checked and recertified by certification bodies and public authorities. In addition our internal auditors monitor our quality in accordance with § 30 of the Act on Medical Devices (MPG). Therefore we improve our overall performance continually.

The product range includes medical sleeping and positioning systems for decubitus therapy, high-quality mattresses and dental applications.

If you would like to find out more about our company and the range of products that we offer, please feel free to contact us at any time – via e-mail, over the phone or even in person.

Further information and news can be found on our homepage funke-medical.de!

We look forward to hearing from you soon!

Andreas Funke
-CEO-

Helmut Wessels
-Executive Board-

Patricia Funke
-Executive Board-



Administration and production facility



Warehouse and Service Centre



Administration and production facility



Environmental responsibility

Energy is a valuable resource that we need for our production.
In order to use energy most efficiently we have taken the following steps:



Since 2017, we generate 250.000 kWh of electricity from solar energy using 4 photovoltaic systems.

The produced amount covers about 65% of our own demand.

This means that approx. 130t less CO₂ will be emitted into the environment every year.



We want to preserve our natural resources with all these energy-saving measures
in order to live up to our environmental responsibility.


**THE WORLD IS A FINE PLACE
AND WORTH FIGHTING FOR.**

Ernest Hemingway



HOMECARE

SYSTEMATIC HOMECARE



Hyper Foam PLUS

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 3
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam PLUS

Decubitus therapy mattress

3 YEARS
GUARANTEE
*²



Size

Height 14 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable



Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 140 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC HOMECARE



Hyper Foam 2

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam 2

Decubitus therapy mattress

5 YEARS GUARANTEE
*²



Size

Height 15 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 150 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC HOMECARE



MAXX 250

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

MAXX 250

Bariatric decubitus therapy mattress (up to 250 kg)

5 YEARS
GUARANTEE
*²



Size

Height 18 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5



Patient weight:
40 - 250 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.


*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



SYSTEMATIC CLINIC CARE



Hyper Foam PLUS^{clinic}

Your retailer:

PRS
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 3
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam PLUS^{clinic}

Decubitus therapy mattress

3 YEARS
GUARANTEE
*²



Size

Height 14 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 140 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



Hyper Foam VISCO_{clinic}

Your retailer:

PRS
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 2
(EPUAP)

Middle risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam VISCO^{clinic}

Decubitus therapy mattress

3 YEARS
GUARANTEE
*²



Size

Height 14 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam +
visco foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 140 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



Hyper Foam 2^{clinic}

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Hyper Foam 2^{clinic}

Decubitus therapy mattress

5 YEARS GUARANTEE
*²



Size

Height 15 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 150 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



MAXX 250
clinic

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

MAXX 250^{clinic}

Bariatric decubitus therapy mattress (up to 250 kg)

5 YEARS
GUARANTEE
*²



Size

Height 18 cm

Length See current price list

Width See current price list



Breathable
open-cell
HR-foam



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Pressure
< 18 mmHg



Bi-elastic cover
to minimise shear
forces



CE - marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5



Patient weight:
40 - 250 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.


*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

Decubitus therapy mattresses with evacuation function

evAQ[®]

- Your **all-in-one** solution -





Decubitus therapy mattresses with evacuation function

- Your **all-in-one** solution -

In case of fire you have only **3 minutes** to escape!



All EVAQ® decubitus therapy systems are equipped with an evacuation function.
Patients can be evacuated by a single person in case of fire.

It takes **less than 60 seconds** to evacuate one person from one floor to the other.

Watch our evacuation film and see for yourself.

Scan me!



FILMLINK

Watch our evacuation film here!




- Your **all-in-one** solution -

Comparison – Furnishing of institutions (e.g. clinics, elderly homes)

Furnishing without EVAQ®



+



Standard care mattress
mostly without
decubitus therapy function

+



Evacuation system
e.g. stretchers,
rescue sheets

In case of evacuation

+




4 – 6 rescuer

Furnishing with EVAQ®



+



=

mattress
+ decubitus therapy system
+ evacuation function

In case of evacuation

+



1 rescuer

Resource saving (money and workforce)!



Decubitus therapy mattresses with evacuation function

All mattresses in this range are equipped with an evacuation function. This means that you can evacuate patients quickly and easily without transferring them on rescue sheets or carry them on stretchers!


Evacuation systems such as stretchers and rescue sheets often require 2 or more persons to use them.

The EVAQ® decubitus therapy mattresses enable you to carry out evacuation gently using only one person, even via stairs!


§§ Quite apart from the fact that every facility operator is responsible for taking measures to ensure rapid evacuation! §§

SAFER → FASTER → GENTLER

- Safety belts with quick-release buckle
- Special sliding material at the bottom
- Stairway evacuation possible
- Firm edges provide stability
- Rescue handles on the cover



SYSTEMATIC CLINIC CARE



evAQ[®] - PRO

- Your **all-in-one** solution -

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy

Up to stage 2
(EPUAP)


Middle risk
(Norton-, Braden-, Waterlow-Scale)

evAQ[®] - PRO

WITH EVACUATION FUNCTION

Decubitus therapy mattress

3 YEARS
GUARANTEE
*²




Size

Height 14 cm

Length See current price list

Width See current price list




*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



evAQ[®] - PRO VISCO

- Your **all-in-one** solution -

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy

Up to stage 2
(EPUAP)


Middle risk
(Norton-, Braden-, Waterlow-Scale)

evAQ® - PRO VISCO

WITH EVACUATION FUNCTION

Decubitus therapy mattress

3 YEARS GUARANTEE
*²




Size

Height 14 cm

Length See current price list

Width See current price list




*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the mattress core (pro rata basis)!

*³ For more details see instructions for use.

SYSTEMATIC CLINIC CARE



evAQ[®] - PRO II

- Your **all-in-one** solution -

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI^{*1} guidelines

Applicable
for pain therapy

Up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

ev**AQ**[®] - PRO II

WITH EVACUATION FUNCTION

Decubitus therapy mattress

5 YEARS
GUARANTEE^{*2}




Size

Height 15 cm

Length See current price list

Width See current price list




^{*1} Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

^{*2} Guarantee applies to functional properties of the mattress core (pro rata basis)!

^{*3} For more details see instructions for use.

SYSTEMATIC CLINIC CARE



evAQ[®] - PRO-XL

- Your **all-in-one** solution -

Your retailer:

PRS 
Pressure Relieving Solutions

Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy

Up to stage 4
(EPUAP)


Very high risk
(Norton-, Braden-, Waterlow-Scale)

evAQ[®] - PRO XL

WITH EVACUATION FUNCTION

Bariatric decubitus therapy mattress (up to 250 kg)

5 YEARS
GUARANTEE^{*2}



Size

Height 18 cm

Length See current price list

Width See current price list



*1 Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*2 Guarantee applies to functional properties of the mattress core (pro rata basis)!


*3 For more details see instructions for use.

DECUBITUS THERAPY SEATING SYSTEMS





Decubitus therapy seat cushions



XAIR®
HYBRID

- 4 For decubitus prophylaxis and therapy up to stage 4 (EPUAP)
- RISK Very high risk (Norton-, Braden-, Waterlow-Scale)




XSEAT®

- 4 For decubitus prophylaxis and therapy up to stage 4 (EPUAP)
- RISK Very high risk (Norton-, Braden-, Waterlow-Scale)




XSEAT®
clinic

- 4 For decubitus prophylaxis and therapy up to stage 4 (EPUAP)
- RISK Very high risk (Norton-, Braden-, Waterlow-Scale)
- Cover is easy to disinfect




XSEAT®
light

- For decubitus prophylaxis
- Cover is easy to disinfect



GEL SEAT®

- 4 For decubitus prophylaxis and therapy up to stage 4 (EPUAP)
- RISK Very high risk (Norton-, Braden-, Waterlow-Scale)



GEL SEAT®
clinic

- 4 For decubitus prophylaxis and therapy up to stage 4 (EPUAP)
- RISK Very high risk (Norton-, Braden-, Waterlow-Scale)
- Cover is easy to disinfect



Funke Medical AG is official sponsoring partner of the German wheelchair basketball national teams

On November 17, 2014 Patricia Funke and Helmut Wessels of the Funke Medical AG and Andreas Joneck, representative for the German Wheelchair Sports Association (Wheelchair Basketball Division) signed a 3-year sponsorship agreement.

basketball teams. Who else could be more authentic?" explained Patricia Funke (head of marketing) and Helmut Wessels (head of sales).

"It was only logical consequence to set up a partnership with the German wheelchair

"We want to create a fruitful cooperation which both sides can benefit from. Therefore our aim is to increase the popularity of this sport and improve the technical possibilities for people having mobility problems." added Helmut Wessels.



German Women's Wheelchair Basketball Team crowned European Champion in 2015.

Wheelchair basketball national player Annabel Breuer is the new prominent face for XSEAT® and GELSEAT® Decubitus Therapy Seating Systems developed by Funke Medical AG.



- 2010: Silver World Championship (Birmingham, GB)
- 2011: Gold European Championship (Nazareth, IL)
- 2012: Gold Summer-Paralympics (London, GB)
- 2013: Silver European Championship (Frankfurt, DE)
- 2014: Silver World Championship (Toronto, CA)
- 2015: Gold European Championship (Worcester, GB)
- 2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
- 2017: Silver European Championship (Tenerife, ES)
- 2018: Bronze World Championship (Hamburg, DE)
- 2019: Bronze European Championship (Rotterdam, NL)


OFFICIAL SPONSOR

German
Wheelchairbasketball
National Teams





DECUBITUS THERAPY SEAT CUSHION



XAIR®
HYBRID

- AIR-CELL SEAT CUSHION -

For active wheelchair users!

Your retailer:

OFFICIAL SPONSOR



2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)

A Breuer
Annabel Breuer
National player
Wheelchair Basketball





Decubitus therapy air-cell seat cushion for active wheelchair users!



Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy

Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Cover:

- Removable
- Washable up to 60°C

Mesh³ fabric guarantees an excellent
microclimate

Waterproof but breathable cover

Individually adjustable
by valve
(one-hand operation)

Unique Geltouch®-technology
with cooling-effect

Hyperelastic three-layered
comfort core:
Geltouch®-top plate

HR-Foam

Non-skid bottom

Size in cm	Item No.	Size in cm	Item No.
Width x Length x Height		Width x Length x Height	
38 x 40 x 8	XAIR3840	45 x 44 x 8	XAIR4544
40 x 40 x 8	XAIR4040	45 x 46 x 8	XAIR4546
40 x 42 x 8	XAIR4042	46 x 51 x 8	XAIR4651
41 x 46 x 8	XAIR4146	48 x 42 x 8	XAIR4842
42 x 40 x 8	XAIR4240	48 x 44 x 8	XAIR4844
42 x 42 x 8	XAIR4242	48 x 46 x 8	XAIR4846
42 x 44 x 8	XAIR4244	50 x 42 x 8	XAIR5042
42 x 46 x 8	XAIR4246	51 x 44 x 8	XAIR5144
42 x 52 x 8	XAIR4252	51 x 46 x 8	XAIR5146
45 x 42 x 8	XAIR4542	52 x 52 x 8	XAIR5252

WATCH NOW
on YouTube!



XAIR® hybrid-Movie



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the cushion core (pro rata basis)!



WATCH NOW
on YouTube!



Handle for easy carrying


Soft protective cover with
hook-and-loop fastener



Lift the protective cover.



Open the valve by pulling the blue
lock and wait approx. 30 seconds,
until the air-cells have inflated
automatically.
(Valve lock can not fall off.)




Close the valve
by pressing.




Put the cushion on the wheelchair or seat.
(The handle must point to the front and the
attached arrows must point upwards.)



Sit on the cushion.



Carefully open the valve again
until the ideal sitting comfort
is achieved.



Close the valve and
the protective cover.



3-chamber seat cushion
provides optimum
pressure relief



Individually adjustable
for ideal seating comfort



No pump or other
tools necessary



Fast and easy
adaptation to external
pressure conditions




No risk of failure due to
self-inflating FOAM-IN-CELL-
TECHNOLOGY



Very light –
70 % less weight than
conventional air-cell seat cushions



DECUBITUS THERAPY SEAT CUSHION



XSEAT®
- SEAT CUSHION -

Your retailer:



2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)

A Breuer
Annabel Breuer
National player
Wheelchair Basketball



Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

XSEAT®


Decubitus therapy seat cushion

3 YEARS GUARANTEE
*²



Size in cm	Item No.
Width x Length x Height	
38 x 40 x 8	X-SEAT3840
40 x 40 x 8	X-SEAT4040
40 x 42 x 8	X-SEAT4042
41 x 46 x 8	X-SEAT4146
42 x 40 x 8	X-SEAT4240
42 x 42 x 8	X-SEAT4242
42 x 44 x 8	X-SEAT4244
42 x 46 x 8	X-SEAT4246
42 x 52 x 8	X-SEAT4252
45 x 42 x 8	X-SEAT4542

Size in cm	Item No.
Width x Length x Height	
45 x 44 x 8	X-SEAT4544
45 x 46 x 8	X-SEAT4546
46 x 51 x 8	X-SEAT4651
48 x 42 x 8	X-SEAT4842
48 x 44 x 8	X-SEAT4844
48 x 46 x 8	X-SEAT4846
50 x 42 x 8	X-SEAT5042
51 x 44 x 8	X-SEAT5144
51 x 46 x 8	X-SEAT5146
52 x 52 x 8	X-SEAT5252




*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the cushion core (pro rata basis)!



DECUBITUS THERAPY SEAT CUSHIONS



XSEAT[®]
clinic

- SEAT CUSHION -

Your retailer:

OFFICIAL SPONSOR



A Breuer
Annabel Breuer
National player
Wheelchair Basketball

2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)



XSEAT[®] clinic

Decubitus therapy seat cushion

3 YEARS
GUARANTEE
*2


Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy

Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Cover is easy to disinfect



Size in cm	Item No.
Width x Length x Height	
38 x 40 x 8	X-SEATCL3840
40 x 40 x 8	X-SEATCL4040
40 x 42 x 8	X-SEATCL4042
41 x 46 x 8	X-SEATCL4146
42 x 40 x 8	X-SEATCL4240
42 x 42 x 8	X-SEATCL4242
42 x 44 x 8	X-SEATCL4244
42 x 46 x 8	X-SEATCL4246
42 x 52 x 8	X-SEATCL4252
45 x 42 x 8	X-SEATCL4542

Size in cm	Item No.
Width x Length x Height	
45 x 44 x 8	X-SEATCL4544
45 x 46 x 8	X-SEATCL4546
46 x 51 x 8	X-SEATCL4651
48 x 42 x 8	X-SEATCL4842
48 x 44 x 8	X-SEATCL4844
48 x 46 x 8	X-SEATCL4846
50 x 42 x 8	X-SEATCL5042
51 x 44 x 8	X-SEATCL5144
51 x 46 x 8	X-SEATCL5146
52 x 52 x 8	X-SEATCL5252




*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the cushion core (pro rata basis)!

*³ For more details see instructions for use.



DECUBITUS THERAPY SEAT CUSHIONS



XSEAT®
light

- SEAT CUSHION -

Your retailer:



2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)

A Breuer
Annabel Breuer
National player
Wheelchair Basketball




Autoclavable according to
RKI*¹ guidelines

For decubitus prophylaxis

Cover is easy to disinfect

XSEAT[®]
light

Decubitus prophylaxis seat cushion



Size in cm Width x Length x Height	Item No.	Size in cm Width x Length x Height	Item No.
38 x 40 x 5	X-SEATL3840	45 x 44 x 5	X-SEATL4544
40 x 40 x 5	X-SEATL4040	45 x 46 x 5	X-SEATL4546
40 x 42 x 5	X-SEATL4042	46 x 51 x 5	X-SEATL4651
41 x 46 x 5	X-SEATL4146	48 x 42 x 5	X-SEATL4842
42 x 40 x 5	X-SEATL4240	48 x 44 x 5	X-SEATL4844
42 x 42 x 5	X-SEATL4242	48 x 46 x 5	X-SEATL4846
42 x 44 x 5	X-SEATL4244	50 x 42 x 5	X-SEATL5042
42 x 46 x 5	X-SEATL4246	51 x 44 x 5	X-SEATL5144
42 x 52 x 5	X-SEATL4252	51 x 46 x 5	X-SEATL5146
45 x 42 x 5	X-SEATL4542	52 x 52 x 5	X-SEATL5252



Breathable
open-cell
HR-foam +
visco foam



Autoclavable



Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Bi-elastic cover
to minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (I+II)
CRIB 5




Patient weight:
20 - 120 kg

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties of the cushion core (pro rata basis)!



DECUBITUS THERAPY SEAT CUSHIONS



WITH
COOLING-
EFFECT!

GEL SEAT[®]
- SEAT CUSHION -

Your retailer:

OFFICIAL SPONSOR



A Breuer
Annabel Breuer
National player
Wheelchair Basketball

2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)



Autoclavable according to
RKI*¹ guidelines

Applicable
for pain therapy


Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)


GEL SEAT®

Decubitus therapy seat cushion

3 YEARS
GUARANTEE
*2




Size in cm	Item No.	Size in cm	Item No.
Width x Length x Height		Width x Length x Height	
38 x 40 x 8	G-SEAT3840	45 x 44 x 8	G-SEAT4544
40 x 40 x 8	G-SEAT4040	45 x 46 x 8	G-SEAT4546
40 x 42 x 8	G-SEAT4042	46 x 51 x 8	G-SEAT4651
41 x 46 x 8	G-SEAT4146	48 x 42 x 8	G-SEAT4842
42 x 40 x 8	G-SEAT4240	48 x 44 x 8	G-SEAT4844
42 x 42 x 8	G-SEAT4242	48 x 46 x 8	G-SEAT4846
42 x 44 x 8	G-SEAT4244	50 x 42 x 8	G-SEAT5042
42 x 46 x 8	G-SEAT4246	51 x 44 x 8	G-SEAT5144
42 x 52 x 8	G-SEAT4252	51 x 46 x 8	G-SEAT5146
45 x 42 x 8	G-SEAT4542	52 x 52 x 8	G-SEAT5252



*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.
*² Guarantee applies to functional properties of the cushion core (pro rata basis)!



DECUBITUS THERAPY SEAT CUSHIONS



GEL SEAT
clinic®

- SEAT CUSHION -

Your retailer:

OFFICIAL SPONSOR



A Breuer
Annabel Breuer
National player
Wheelchair Basketball

2010: Silver World Championship (Birmingham, GB)
2011: Gold European Championship (Nazareth, IL)
2012: Gold Summer-Paralympics (London, GB)
2013: Silver European Championship (Frankfurt, DE)
2014: Silver World Championship (Toronto, CA)
2015: Gold European Championship (Worcester, GB)
2016: Silver Summer-Paralympics (Rio de Janeiro, BR)
2017: Silver European Championship (Tenerife, ES)
2018: Bronze World Championship (Hamburg, DE)
2019: Bronze European Championship (Rotterdam, NL)



GEL SEAT
clinic®

Decubitus therapy seat cushion

3 YEARS
GUARANTEE
*2

Autoclavable according to
RKI*1 guidelines

Applicable
for pain therapy

Up to stage 4
(EPUAP)

Very high risk
(Norton-, Braden-, Waterlow-Scale)

Cover is easy to disinfect

PU Cover:

- Removable
- Washable up to 95°C
- Disinfection by wiping
- Disinfectant-proof*3

Ultrasonically welded
seams

Waterproof but breathable cover

Anatomic bodyshape
ensures a good posture
in a seated position


Zipper cover

Non-skid bottom

Hyperelastic two-layered
comfort core:
Geltouch®-top plate

Unique Geltouch®-technology
with cooling-effect

Size in cm	Item No.	Size in cm	Item No.
Width x Length x Height		Width x Length x Height	
38 x 40 x 8	G-SEATCL3840	45 x 44 x 8	G-SEATCL4544
40 x 40 x 8	G-SEATCL4040	45 x 46 x 8	G-SEATCL4546
40 x 42 x 8	G-SEATCL4042	46 x 51 x 8	G-SEATCL4651
41 x 46 x 8	G-SEATCL4146	48 x 42 x 8	G-SEATCL4842
42 x 40 x 8	G-SEATCL4240	48 x 44 x 8	G-SEATCL4844
42 x 42 x 8	G-SEATCL4242	48 x 46 x 8	G-SEATCL4846
42 x 44 x 8	G-SEATCL4244	50 x 42 x 8	G-SEATCL5042
42 x 46 x 8	G-SEATCL4246	51 x 44 x 8	G-SEATCL5144
42 x 52 x 8	G-SEATCL4252	51 x 46 x 8	G-SEATCL5146
45 x 42 x 8	G-SEATCL4542	52 x 52 x 8	G-SEATCL5252



*1 Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*2 Guarantee applies to functional properties of the cushion core (pro rata basis)!


*3 For more details see instructions for use.

POSITIONING AIDS





SYSTEMATIC DECUBITUS THERAPY



Your retailer:

POSIMED®
- POSITIONING AIDS -

Autoclavable according to
RKI^{*1} guidelines

Applicable
for pain therapy

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED[®]

Universal Cushion

5 YEARS
GUARANTEE^{*2}

Cushion:

- Disinfection by wiping
- Autoclavable

Biocompatible

Waterproof but breathable
incontinence cover ensures an
excellent microclimate

Cushion filling ensures:

- Rapid adaptation to the contours of the body
- Optimised pressure relief

POSIMED[®] positioning aids are available in many different shapes and sizes.

Item description	Size in cm Length x Width	Item No.
Universal Cushion small	37 x 26	PO-UNIK3726
Universal Cushion large	56 x 40	PO-UNIG5640
Heel Rail	39 x 21	PO-FERS3921
Forearm Splint	39 x 21	PO-UNTE3921
Abduction Pillow	60 x 35	PO-SPRE6035

Item description	Size in cm Length x Width	Item No.
Bolster	60 x 21	PO-ROLL6021
30° Positioning Cushion small	170 x 74	PO-LAGE0165
30° Positioning Cushion large	182 x 54	PO-LAGE1165
Hip Abduction Pillow	35 x 28	PO-HÜFT3528
Semi-Circle	137 x 80	PO-HALB1378



Autoclavable




Cover material virus-
and bacteria-proof
accord. ASTM
1670 + 1671



Waterproof but
permeable to
water vapour



Optimised
pressure relief



Bi-elastic cover
to minimise shear
forces



CE marking
MDR 2017/745



Maintenance-free



Flame retardant
DIN EN 597 (H-II)
CRIB 5



Do not wash!



Easy to clean by
wipe disinfection



POSIMED® Product Range

Universal Cushion



e.g. for multi-functional positioning aid for various purposes

Size in cm	Item No.
37/26	PO-UNIK3726
56/40	PO-UNIG5640

Heel Rail



e.g. for positioning the heel

Size in cm	Item No.
39/21/8	PO-FERS3921

Forearm Splint



e.g. for positioning the forearm

Size in cm	Item No.
39/21/8	PO-UNTE3921

Abduction Pillow



e.g. for positioning the knee, ankle and epicondyle

Size in cm	Item No.
60/35	PO-SPRE6035

Bolster



e.g. for supporting and positioning the knee, neck and other parts of the body

Size in cm	Item No.
60/21	PO-ROLL6021

Hip Abduction Pillow



e.g. for positioning the hip
→ reduces the risk of dislocating hip prosthesis

Size in cm	Item No.
35/28/14	PO-HÜFT3528

30° Positioning Cushion



e.g. for support in the stable 30° position, with integrated pillow

Size in cm	Item No.
170/74	PO-LAGE0165
182/54	PO-LAGE1165

Semi-Circle



e.g. for upper body support

Size in cm	Item No.
137/80	PO-HALB1378



POSIMED®

Positioning aids

POSIMED® cushions are positioning aids in many different shapes and sizes for pressure sore prevention and treatment support.


With the various positioning aids the patient can be moved marginally and the well-balanced pressure distribution enables the patient to lie in a relaxed and pressure-relieving position.

The aim is to keep the duration and the intensity of pressure to a minimum.


Our positioning system is based on a pressure-relieving and malleable ball-shaped filling.

Positioning examples


30° lateral position




Heel protection



Semi-sitting position



135° lateral position






Systematic Decubitus Therapy



 **POSIMED®**
HEEL PROTECT®
- Heel protector boot -

Your retailer:



Autoclavable according to
RKI*¹ guidelines

Suitable for the treatment
of diabetic foot ulcers

For decubitus prophylaxis
and therapy

Cover is easy to disinfect

POSIMED® **HEEL PROTECT®**

Heel protector boot

1 YEAR
GUARANTEE*²

Soft fastener

Soft fastener for
drop foot prophylaxis



Mesh³ fabric for
maximum comfort
and climate control

Cover is easy to disinfect:

- Washable up to 95°C
- Suitable for tumble dryer
- Autoclavable
- Disinfection by wiping

ROTATION-STOP



Unique
GELTOUCH®-technology
with cooling-effect

Size

Item No.

ONE SIZE FITS ALL

PO-HEELPR




Breathable
open - cell
Geltouch®-foam +
HR-foam



Autoclavable



Optimised
pressure relief



Bi-elastic cover
to minimise shear
forces



CE marking
MDR 2017/745



Flame retardant
DIN EN 597 (I+II)
CRIB 5



Designed to fit
a person of
any size

*¹ Robert-Koch-Institute (RKI) is a part of the Federal Ministry of Health in Germany. It is responsible for disease prevention and control.

*² Guarantee applies to functional properties (pro rata basis)!



www.funke-medical.de