

Juzo Compression Wrap

Size charts

Juzo Compression Wrap · Foot

Art. 6000 LF/LFL	Standard sizes · Circumferences in cm				
	2/S	3/M	4/L	5/XL	
	15.0–23.0	20.0–28.0	25.0–33.0	30.0–40.0	
	19.0–24.0	24.0–29.0	29.0–34.0	34.0–40.0	
	normal	16	18	19.5	19.5
	long	19	21	22.5	23

* 6 cm proximal to the little toe metatarsophalangeal joint
 ** The circumference can be adjusted using the extended tab at “cB” if required.

Juzo Compression Wrap · Calf

Art. 6000 LC/LLC/LCM/LLCM	Standard sizes · Circumferences in cm						
	2/S	3/M	4/L	5/XL	6/XXL	7/XXXL	
 pD AD normal < 40.0 long > 40.0	cC	29.0–35.0	34.0–40.0	37.0–44.0	41.0–48.0	46.0–53.0	50.0–60.0
	max	35.0–41.0	40.0–46.0	44.0–51.0	48.0–55.0	53.0–60.0	60.0–70.0
	cB	18.0–21.0	21.0–24.0	24.0–27.0	27.0–31.0	31.0–35.0	36.0–47.0

Juzo Compression Wrap · Knee

Art. 6000 LK	Standard sizes · Circumferences in cm						
	2/S	3/M	4/L	5/XL	6/XXL	7/XXXL	
	cE	30.0–36.0	36.0–42.0	42.0–48.0	48.0–55.0	55.0–62.0	62.0–70.0



Juzo Compression Wrap · Thigh

Art. 6000 LT/LTg	Standard sizes · Circumferences in cm					
	2/S	3/M	4/L	5/XL	6/XXL	
 pFG normal < 25.0 long > 30.0	cG	43.0–70.0	54.0–81.0	65.0–92.0	75.0–102.0	85.0–112.0
	cF	35.0–55.0	48.0–70.0	60.0–82.0	70.0–92.0	80.0–102.0




Juzo Compression Wrap

Size charts

Juzo Compression Wrap · Hand

Art. 6000 UHR/UHL	Standard sizes · Circumferences in cm				
	1/XS	2/S	3/M	4/L	5/XL
 cA	15.5–17.5	17.5–19.5	19.5–22.0	22.0–24.5	24.5–28.0
 cC	13.0–17.0	15.0–19.0	17.0–21.5	20.0–24.0	21.0–25.0

Juzo Compression Wrap · Arm

Art. 6000 UA/UAG Lengths in cm	Standard sizes · Circumferences in cm				
	1/XS	2/S	3/M	4/L	5/XL
 cG	27.0–36.0	30.0–39.0	33.0–44.0	35.0–49.0	40.0–56.0
 cE	23.0–32.0	24.0–34.0	26.0–36.0	29.0–42.0	31.0–46.0
 cC	16.0–19.0	17.0–20.0	18.0–22.0	20.0–26.0	21.0–29.0

PCG
normal < 43,0
long > 43,0

Juzo Liner

Art. 6082

Liner stocking for medical adaptive compression systems

- protects sensitive skin
- ideal for securing wound dressings when the product is put on or taken off
- perfect for use with medical adaptive compression systems

Characteristics

Complementary products



Application area

Phlebology
Lymphology

Features

+ ideal for fixation

Treatment area

+ lower limb

Size chart

Art. 6082 Lengths in cm	Standard sizes · Circumferences in cm	
	I/XS	II/S
PG AG normal 80.0	50.0–75.0	76.0–99.0
PD AD normal 45.0	29.0–48.0	48.0–69.0

Lower limb

Colour

Black



	Form	Model	Unit	Standard sizes
		Can only be used as underliner stocking for medical adaptive compression systems		6082 €
	AD	Below-knee stockings	Pair	Size 1–2 36.40
	AG	Thigh stockings	Pair	Size 1–2 48.60

Product examples



AD · Black
Juzo Liner



AG · Black
Juzo Liner