

Adjusted by:

E-mail:

Adjusted for:

Tel:

Date:

Total weight (g): **0**

The weight changes based on the choices made below.

S3 Long has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". S3 has a seat height of 45 cm in the front and 43 at the rear, without cushion. To adjust the seat height you can choose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options. **Max userweight: 100 kg for seat width 33-42 and 150kg for seat width 45-50 cm.**

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. **BB** = seatwidth, **HH**= backrest height

Original design

Alternative options

Other options

[info](#)



1) Chassis-colour

Click on the link to see up-to-date colours.

Fill in art.no. in the square

2) Seat width (BB)

[art.no.](#)

Clear

- 33 cm 36 cm 39 cm 42 cm 45 cm 50 cm

Art.no include: chassis, seat, calf strap, backrest lock

3) Seat depth

[art.no.](#)

Clear

45 cm

42,5 cm

4) Backrest

[Standard](#)
523BBHH VER2

Clear

Height (HH)

- 20 cm 35 cm
 25 cm 40 cm
 30 cm 45 cm

Clear

[Tapered](#)

525BBHH VER2

[Curved back](#)

524BBHH VER2

[Backrest, extended range](#)

[Frame for fixed backrest systems](#)

art. _____

5) Rear axle

[Carbon standard](#)
33023 BB

Clear

[Carbon +15 mm](#)

33023BC*

* seat width **BB** +1

Use this rear axle when choosing arm support at Accessories and/or increased seat height

[Other](#) rear axle or seat height

[Increase the seat height](#) +0 mm

6) Castor and fork

[120 mm](#)

Clear

[90 mm](#)

[150 mm](#)

7) Footrest

[Footrest Titanium](#)
36001BB

Clear

[Extended, reinforced](#)

36111BB

[Foldable footplates](#)

[art.no.](#)

[Heel strap for footplates](#)

[art.no.](#)

[Plastic plate for footrest](#)

36010BB

[Other footrest](#)

8) Brake

[S3 standard](#)
4630000

Clear

[One arm brake, right](#)

[art.no.](#)

[One arm brake, left](#)

[art.no.](#)

9) Rear wheel - without pushrim

[24" standard](#)
[art.no.](#)

Clear

[Rear wheel sizes](#)

- 24" [Spinergy X](#) black spokes [art.no.](#)
 24" [Spinergy X](#) white spokes
 25" [Spinergy X](#) black spokes
 25" [Spinergy X](#) white spokes
 26" [Spinergy X](#) black spokes
 26" [Spinergy X](#) white spokes

[Other rear wheel](#)

Tire - choose tire

[Schwalbe Right Run](#) 24",25",26"

Clear

[Rugged](#) 24"

Schwalbe [Marathon Plus](#) 24",25",26"

Schwalbe [One](#) 24",25"

[PU](#) 24" (massive)

[reTyre](#)-tyre skins (does not fit PU tyres)

10) Pushrim

[Titanium](#)
190DD11 (x2)

Clear

DD = Wheel diameter

[Alu X](#) 24"

1904404 (x2)

[Alu X](#) 25"

1904405 (x2)

[Alu X](#) 26"

1904406 (x2)

[Tetragrip](#)

19052DDBT x2

[Paragrip](#)

19051DDBT x2

DD = Wheel diameter

[Maxgrip](#) 24" (contains latex)
1903060 x2

[Maxgrip](#) 22" (contains latex)
1903080 x2



Accessories



Side guard

- [Side guard S3](#) [Video](#)
[art.no.](#)
- Clear



Side guard Carbon

- [Carbon](#)
[art.no.](#)
 - [Carbon, small cover](#)
[art.no.](#)
 - [Carbon, large cover *](#)
[art.no.](#)
- * Only with one arm brake



Arm support

- [Arm support S3](#) [Video](#)
[art.no.](#)



Push handles H/S

- [Height adjustable S3](#) [Video](#)
4252060 VER2
- [Height adjustable S3+10 cm](#)
4252063 VER2
- Max. userweight 100kg
- Clear



Push handle Foldable

- [Foldable S3](#) [Video](#)
4245000



Push bar

- [Push bar S3](#)
42400BB VER2
- Max. userweight 100kg



Anti tip

- [Anti tip](#) std seat height
4353000 [Video](#)
- [Other seat height](#)



Cushion

- Seat depth 35cm Seat depth 40cm
- [2,5 cm](#) [2,5 cm](#)
 - 572BB35 572BB40
 - [5 cm](#) [5 cm](#)
 - 575BB35 575BB40
 - Clear



Spoke protectors

- [24" logo](#)
4500000
- [24" no logo](#)
4500011
- Clear




Other

- [Tube protection](#) _____ p
9900251 (x3/wheelchair)
- [Crutch holder S3](#) _____ p
9900301 (x1/crutch) [Video](#)
- [Backrest wedge](#) _____ p
4710012

Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)

- [Freewheel](#) Height 6-14 cm art.no. 2600102



Wedge cushion

- [Wedge cushion](#)
57710BB

Set-up adjustments



Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=60 mm


A= _____ mm



Footrest height [Footbow- Footplates](#)

Note measure (A) and (B) for desired footrest.
[Bambino foot plate with Feewheel - info](#)

A= _____ mm B= _____ mm



Transportation marking

- YES do mark the wheelchair with attachment points for transport in disabled service vehicle

Art.no.9000152

NOTE! See [manual](#) for more information.



Extra