

Adjusted by:

E-mail:

Adjusted for:

Tel.:

Date:

Total weight (g): **0**

The weight changes based on the choices made below.

U3 Y-front low has a adjustable rear axle in carbon fibre. To find a balance that suits you the rear axle can be adjusted into different positions. Note desired measure at "Settings". U3 Y-front low has a seat height of 44,5 cm in the front and 41,5 at the rear, without cushion. To adjust the seat height you can chose other rearwheel sizes or a raising Kit with suitable castor fork from the List of options.

Max userweight: 100 kg for seat width 33-42.

Save the guide to the "desktop" and open it with [Adobe Reader](#), for best function Design your wheelchair by clicking in the circles in this guide. Blue, underlined text is clickable to [show info](#) or image. BB = seatwidth, HH= backrest height

Original design

Alternative options

Other options

[info](#)

1) Chassis-colour [Click on the link to see up-to-date colours](#)

Choose colour



2) Seat width (BB)

[art.no.](#)

Clear

- 33 cm 36 cm 39 cm 42 cm 45 cm 50 cm

Art.no include: chassis, seat, calf strap, backrest lock



3) Seat depth

[art.no.](#)

Clear

35 cm

30 cm

37,5 cm

Only alternative seat at seat width 50 cm



4) Backrest

[Standard](#)
523BBHH VER2

Clear

Height (HH)

- 20 cm 35 cm
 25 cm 40 cm
 30 cm 45 cm

Clear

[Tapered](#)
525BBHH VER2

[Curved back](#)
524BBHH VER2

[Backrest, extended range](#)

[Frame fo fixed backrest systems](#)

art. _____



5) Rear axle

[Carbon standard](#)
33023BB

Clear

[Carbon +15 mm](#)
33023BC*

* seat width BB +1

Use this rear axle when chosing, arm support at Accessories and/or increased seat height

[Other](#) rear axle or seat height

[Increase the seat height](#) +0 mm



6) Castor and fork

[90 mm](#)
[art.no.](#)

Clear

[120 mm](#)
[art.no.](#)



7) Footrest [info about footrest height](#)

[Footrest Titanium](#)
36915BB

Clear

[Extended, reinforced](#)
36006BB, needed for Freewheel

[Plastic plate for footrest](#)
36010BB

[Other footrest](#)



8) Brake

[Brake S3](#)
4630000

Clear

[One arm brake, right](#)
[art.no.](#)

[One arm brake, left](#)
[art.no.](#)



9) Rear wheel - without pushrim

24" [standard](#)
[art.no.](#)

Clear

- [Rear wheel sizes](#) [art.no.](#)
- 24" [Spinergy X](#) black spokes
 - 24" [Spinergy X](#) white spokes
 - 25" [Spinergy X](#) black spokes
 - 25" [Spinergy X](#) white spokes
 - 26" [Spinergy X](#) black spokes
 - 26" [Spinergy X](#) white spokes

[Annat drivhjul](#)



Tire - choose tire

[Schwalbe Right Run](#) 24",25",26"

Clear

- [Rugged 24"](#)
- Schwalbe [Marathon Plus](#) 24",25",26"
- Schwalbe [One](#) 24",25"

[PU 24"](#) (massive)

[reTyre](#)-tyre skins (does not fit PU tyres)



10) Pushrim

[Titanium](#)
190DD11 (x2)

Clear

DD = Wheel diameter

- [Alu X 24"](#)
1904404 (x2)
- [Alu X 25"](#)
1904405 (x2)
- [Alu X 26"](#)
1904406 (x2)
- [Tetragrip](#)
19052DDBT x2
- [Paragrip](#)
19051DDBT x2

DD = Wheel diameter

[Maxgrip 24"](#) (contains latex)
1903060 x2


Accessories

Side guard *

[Side guard S3](#) [Video](#)
[art.no.](#)

* Use sideguard 26" for 24" rear wheels due to the low chassis.

Clear



Side guard Carbon

- [Carbon](#)
[art.no.](#)
- [Carbon, small cover](#)
[art.no.](#)
- [Carbon, large cover *](#)
[art.no.](#)

Use sideguard 26" for 24" rear wheels because of the low chassis.

* Only with one arm brake



Arm support *

[Arm support S3 *](#) [Video](#)
[art.no.](#)

* Not possible to combine with 25" and 26" rear wheels
Use Arm support 26" for 24" rear wheels due to the low chassis

Push handles H/S

[Height adjustable S3](#) [Video](#)
4252060 VER2

[Height adjustable S3+10 cm](#)
4252063 VER2
Max. userweight 100kg

Clear



Push handle Foldable

[Foldable S3](#) [Video](#)
4245000



Push bar

[Push bar S3](#)
42400BB VER2

Max. userweight 100kg

Anti tip low chassis

[Anti tip](#) std seat height
4355000 [Video](#)

[Other seat height](#)



Cushion

Seat depth 35cm Seat depth 40cm

[2.5 cm](#) [2.5 cm](#)

572BB35 572BB40

[5 cm](#) [5 cm](#)

575BB35 575BB40

Clear



Spoke protectors

[24" logo](#)

4500000

[24" no logo](#)

4500011

Clear

Other

[Tube protection](#) _____ p
9900251 (x3/wheelchair)

[Crutch holder S3](#) _____ p
9900301 (x1/crutch)
[Video](#)


[Backrest wedge](#) _____ p
4710012

Chose art.no. for Freewheel according to the **footrest height** of the wheelchair Read more about Freewheel on the [homepage](#)

[Freewheel](#) Height 6-14 cm art.no. 2600102

Wedge cushion

[Wedge cushion](#)
57710BB



Set-up adjustments

Back rest angle

Note measure (A) for desired back rest angle.
14-28 mm. Standard=20 mm

A= _____ mm



Rear axle adjustment

Note measure (A) for desired balance.
Standard=60 mm

A= _____ mm



Footrest height [info about footrest height](#)

Note measure (A).

A= _____ mm



Transportation marking

YES do mark the wheelchair with attachment points for transport in disabled service vehicle
Art.no.9000152

NOTE! See [manual](#) for more information.



Extra