

SmartDrive® MX2+ with PushTracker

SmartDrive is a lightweight, power assist for manual wheelchairs that allows you to go up hills and ramps, over thick carpet, long distances, and essentially, anywhere where it is hard to push.

It is simple and intuitive to operate, with the help of the PushTracker BLUETOOTH® wristband.

- Double tap to start. You can start from a dead stop or activate SmartDrive in mid-roll
- A single tap at coasting speed maintains the speed you are at – no need to push.
- When you are ready to stop, double tap again and then brake with your hands.

Simple!

Increasing your independence

SmartDrive greatly reduces the amount of times you have to push in a day. This helps to reduce fatigue and pain in shoulders and arms, which means you can get out and about more on your own.

"I can do more and not limit my life, because I don't have to worry as much about whether I'll have the strength and energy to push through the whole day.

I can carry things, like a bag of shopping, and not have to think about how I'm going to propel the chair. SmartDrive gives you the energy to do more – which means you get to experience more in life."

Chelsea Smith, Nashville, USA





Keeping you healthy

Research shows that many manual wheelchair users experience upper extremity pain. SmartDrive helps you stay healthy by protecting your shoulders from wear and tear caused by long-time, repetitive propulsion and reduces the risk for RSI (repetitive strain injury) on your elbows and wrists.

Because it is light and compact, it can be easily mounted or removed. No need for unnecessary transfers or lifting. All of these things save your energy for the important activities of your daily life.





Specifications

Range (level/no resistance conditions)	19.8 km
Maximum speed (Level surface)	6 km/h (default), adjustable
Maximum speed (6° incline):	6 km/h (default), adjustable
Maximum acceleration:	Adjustable up to 0.6m/s2
	(0,2 default)
Motor power:	250W
Motor voltage:	36V
User weight range:	14kg to 150 kg
Operating temperature:	-25°C to 50°C
Sound level	60.7 dB (A)
Total weight:	6.1 kg
Expected lifetime/service life:	5 years
Expected lifetime/service life: Contains BLUETOOTH LE: FCC ID	5 years QOQBT113
Contains BLUETOOTH LE: FCC ID	QOQBT113
Contains BLUETOOTH LE: FCC ID	QOQBT113 5123A-BGTBLE113
Contains BLUETOOTH LE: FCC ID IC ID Overall length:	QOQBT113 5123A-BGTBLE113 389 mm
Contains BLUETOOTH LE: FCC ID IC ID Overall length: Overall width:	QOQBT113 5123A-BGTBLE113 389 mm 141 mm
Contains BLUETOOTH LE: FCC ID IC ID Overall length: Overall width: Overall height:	QOQBT113 5123A-BGTBLE113 389 mm 141 mm 242 mm
Contains BLUETOOTH LE: FCC ID IC ID Overall length: Overall width: Overall height: Omni-wheel diameter:	QOQBT113 5123A-BGTBLE113 389 mm 141 mm 242 mm
Contains BLUETOOTH LE: FCC ID IC ID Overall length: Overall width: Overall height: Omni-wheel diameter: Omni-wheel width:	QOQBT113 5123A-BGTBLE113 389 mm 141 mm 242 mm 193 mm 71 mm

NOTE: The stated values are from testing with a 150kg (331lbs) user. The range, speed, etc. can vary depending on user weight and driving conditions.

Wheelchair compatibility

The MX2+ is compatible with any manufacturer's rigid or folding frame manual wheelchair for users from 14 to 150 kg. It can be mounted to any class of wheelchair including standing chairs and tilt in space.

Wheelchair frame type:	Rigid
	Folding
	Tilt-in-Space
	One-Arm Drive
	Standing
Drive wheel diameter:	501, 540, 559 or 590 mm
Folding frame attachment distance:	210 - 533 mm



