

Valid from 090401

Last change 120824

Article number 803380

### **SEAT SYSTEMS**

For the REAL chairs











#### **ERGOMEDIC - ADVANCED ERGONOMICS**

With various backs, seats, castors, lifting columns and accessories, it is easy to put a functional, comfortable chair together. To further increase sitting comfort, Mercado Medic has collaborated with registered physiotherapist, Bengt Engström, to develop the ErgoMedic seat system for the REAL 9000. ErgoMedic consists of three interacting parts: Seat, Back and the Armrest.

#### **Backrest**

The ERGOMEDIC back rest and seat work together to provide the most comfortable, and effective sitting position. The shape and angle of the high back has been improved to create more space for the pelvis and its softer parts. The lower edge of the back panel now has a gentle 'U' shape and is softer. This, combined with the shape of the seat, provides better sitting depth and a more comfortable position for the lower back and pelvis.

#### Seat

For optimal sitting comfort, the ERGOMEDIC® seats have the right form and softness to ensure optimal stability and comfort, higher at the thigh than the buttocks. Its form distributes pressure evenly and effectively, and prevents the user from sliding forward when the seat is angled. The back edge of the seat is set low to increase the sitting depth.

This means that the seat bone nodules can move backwards when the users stretch out forwards. When the pelvis is positioned at the very back of the seat, and well supported by the backrest, the users legs are also stabilised.

#### Coxit seat


Mercado's Coxit seat has the same high ergonomic quality as the ERGOMEDIC seat. The Coxit seat is divided with individually adjustable flaps. It is used primarily for the alleviation and rehabilitation of orthopaedic injuries to hips or thighs.

#### **Armrest**

The armrest is an extremely important part of the ERGO-MEDIC system, so we have also developed these to be more ergonomic and comfortable. The edges have been bevelled, both on the inside and the outside to minimise the pressure against the underarms during prolonged contact.

The armrest panels are hollowed out at the back in order to stabilise the elbows. At the front they have been designed to provide a good grip when standing up and transferring pressure.

The armrests are covered with soft imitation leather that is more pleasant to the touch and is less slippery if the armrest becomes damp.



| Bac   | krest ErgoMedic, width x depth   | Article No. |  |  |  |
|---|--|-------------|--|--|--|
| 1   | Backrest low ErgoMedic 27x22 cm  | BR107046    |  |  |  |
| 2   | Backrest low ErgoMedic 30x26 cm  | BR106546    |  |  |  |
| 3   | Backrest low ErgoMedic 33x26 cm  | BR106046    |  |  |  |
| 4   | Backrest low ErgoMedic 38x26 cm  | BR105546    |  |  |  |
| 5   | Backrest low ErgoMedic 27x32 cm  | BR104546    |  |  |  |
| 6   | Backrest low ErgoMedic 30x34 cm  | BR103546    |  |  |  |
| 7   | Backrest low ErgoMedic 33x36 cm  | BR102546    |  |  |  |
| 8   | Backrest low ErgoMedic 40x34 cm  | BR103046    |  |  |  |
| 9   | Backrest low ErgoMedic 42x43 cm  | BR101546    |  |  |  |
| Sea   | it ErgoMedic, width x depth  | Article No. |  |  |  |
| 10  | Seat ErgoMedic 29x32 cm  | ST105546    |  |  |  |
| 10  | Seat ErgoMedic 33x35 cm  | ST105046    |  |  |  |
| 10  | Seat ErgoMedic 36x39 cm  | ST104546    |  |  |  |
| 10  | Seat ErgoMedic 36x43 cm  | ST104046    |  |  |  |
| 10  | Seat ErgoMedic 40x39 cm  | ST103546    |  |  |  |
| 10  | Seat ErgoMedic 40x43 cm  | ST103046    |  |  |  |
| 10  | Seat ErgoMedic 40x48 cm  | ST102546    |  |  |  |
| 10  | Seat ErgoMedic 44x43 cm  | ST102046    |  |  |  |
| 10  | Seat ErgoMedic 44x48 cm  | ST101546    |  |  |  |
| 10  | Seat ErgoMedic 48x48 cm  | ST101046    |  |  |  |
| Seat Coxit ErgoMedic, width x depth Article No. |  |             |  |  |  |
| 11  | • •  | ST133046    |  |  |  |
| 11  | Seat Coxit Int ErgoMedic 40x43 cm<br>Seat Coxit Int ErgoMedic 40x48 cm | ST132546    |  |  |  |
| 11  | Seat Coxit Int ErgoMedic 40x46 cm                                      | ST132046    |  |  |  |
| 11  | Seat Coxit Int ErgoMedic 44x48 cm                                      | ST132046    |  |  |  |
| 11  | Seat Coxit Int ErgoMedic 44x48 cm                                      | ST131046    |  |  |  |
| 1.1   | Jeal John Int Ligowieur 40146 Cm                                       | 31131040    |  |  |  |

EXPLANATION OF THE STRUCTURE OF THE ARTICLE NUMBER

Article No: ST101546

"ErgoMedic seat, plush grey 44x48 cm"

Digits 1-2 = Type of upholstered article Digits 3-4 = Seat system
Digits 5-6 = Size
Digit 7 = Material
Digit 8 = Colour



# ERGOMEDIC® PLUS – ERGONOMIC DESIGN THAT IS OUT OF THE ORDINARY

#### **Backrest**

The ERGOMEDIC® PLUS back is available in several variants. The higher backrest provides more shape around the pelvis and thoracic spine, whereas the lower backrest has the same shape, but supports the pelvis and the lumbar spine.

The lower part of the high backrest is shaped so that you can sit with pelvic support even when tilted backwards. At the top in the centre, between the shoulder blades, a soft extra support has been integrated. This to enable the user to sit upright and stable against the backrest, whilst at the same time leaving enough room for free movement of the shoulders.

The high backrest is also available with extra side wedges. The low backrest has a similar shape to the high one, but is designed for supporting the pelvis and lumbar spine.

#### Seat

The biggest requirement of the seat is to be stable and comfortable. ERGOMEDIC® PLUS seats have the right form and softness to ensure optimal stability and comfort, higher under the thigh than the buttocks. Its form distributes pressure evenly and effectively, and prevents the user from sliding forward when the seat is angled. When the user sits at the back of the seat, so that the pelvis is firmly placed against the lower parts of the back support, their legs are also

stabilized. The higher rounded shape on either side of the thighs means that both the pelvis and the legs are stabilized.

#### Armrest

The user should be able to switch between active, and relaxed positions. The ERGOMEDIC® PLUS armrests are shaped so the arms are supported from the sides, which means that the upper body is more relaxed. They are also raised beneath the wrists to provide extra relief.

At the front, the armrest has a rounded convex shape, which makes it comfortable to rest the hand on , while providing a good grip when you want to get in or out of the chair

#### **Neckrest**

The neckrest provides support primarily against the rear, lower part of the head. The lower part of the neck rest is built up with soft padding on

each side, providing support for the lower skull and neck. The upper part is flatter, which leaves room for the user to turn their head. The shape and properties of the neck rest means that the head is well stabilised, even when the user is sitting upright.



| Bac | krest ErgoMedic PLUS, width x depth    | Article No. | Sea     | t ErgoMedic PLUS, width x depth             | Article No. |
|-----|--|-------------|---------|---|-------------|
| 1   | Backrest ErgoMedic Plus 35x20 cm       | BR208046    | 9       | Seat ErgoMedic Plus 40x43 cm                | ST203046    |
| 1   | Backrest ErgoMedic Plus 40x20 cm       | BR207546    | 9       | Seat ErgoMedic Plus 40x48 cm                | ST202546    |
| 2   | Backrest ErgoMedic Plus 30x34 cm       | BR203546    | 9       | Seat ErgoMedic Plus 44x43 cm                | ST202046    |
| 3   | Backrest ErgoMedic Plus 35x43 cm       | BR202046    | 9       | Seat ErgoMedic Plus 44x48 cm                | ST201546    |
| 4   | Backrest ErgoMedic Plus 42x43 cm       | BR201546    | 9       | Seat ErgoMedic Plus 48x48 cm                | ST201046    |
| 4   | Backrest ErgoMedic Plus 46x43 cm       | BR201046    | •       | Geat E. gemeale i iae ien ie em             | 0.20.0.0    |
| 5   | Backrest EM Plus w wedge 30x34 cm      | BR253546    | Sea     | t Coxit ErgoMedic PLUS, width x depth       | Article No. |
| 6   | Backrest EM Plus w wedge 35x43 cm      | BR252046    | 10      | Seat Coxit Int EM Plus 40x39 cm             | ST233546    |
| 6   | Backrest EM Plus w wedge 42x43 cm      | BR251546    | 10      | Seat Coxit Int EM Plus 40x43 cm             | ST233046    |
| 6   | Backrest EM Plus w wedge 46x43 cm      | BR251046    | 10      | Seat Coxit Int EM Plus 40x48 cm             | ST232546    |
| 7   | Backrest EM Plus reinf wedge 35x43 cm  | BR352046    | 10      | Seat Coxit Int EM Plus 44x43 cm             | ST232046    |
| 7   | Backrest EM Plus reinf wedge 42x43 cm  | BR351546    | 10      | Seat Coxit Int EM Plus 44x48 cm             | ST231546    |
| 7   | Backrest EM Plus reinf wedge 46x43 cm  | BR351046    | 10      | Seat Coxit Int EM Plus 48x48 cm             | ST231046    |
| 8   | Backrest EM Plus Comfort 35x43 cm      | BR212046    |         | Cour Coxie in Livinido Toxio on             | 0.20.0.0    |
| 8   | Backrest EM Plus Comfort 42x43 cm      | BR211546    | Nec     | krest, width x height                       | Article No. |
| 8   | Backrest EM Plus Comfort 46x43 cm      | BR211046    | 11      | Nekrest cushion EM Plus 30x23 cm            | NS2046      |
| 8   | Back EM Plus Comfort w wedge 35x43     | BR262046    | 12      | Nekrest cushion EM Plus 30x8 cm             | NS4046      |
| 8   | Back EM Plus Comfort w wedge 42x43     | BR261546    | '-      | Tronical addition Livin tae coxe on         | 110 10 10   |
| 8   | Back EM Plus Comfort w wedge 46x43     | BR261046    | Armrest |   | Article No. |
| 8   | Back EM Plus Comfort reinf wedge 35x43 | BR362046    | 13      | Armrest ErgoMedic Plus Right 42 cm          | AR205138    |
| 8   | Back EM Plus Comfort reinf wedge 42x43 |             | 13      | Armrest ErgoMedic Plus Left 42 cm           | AR205238    |
| 8   | Back EM Plus Comfort reinf wedge 46x43 | BR361046    | 13      | Armrest ErgoMedic Plus Right 30 cm          | AR202138    |
| Ü   | Back EM Flas Comfort form weage 40x40  | D1(001040   | 13      | Armrest ErgoMedic Plus Left 30 cm           | AR202238    |
| Sea | t ErgoMedic PLUS, width x depth        | Article No. | 10      | Anniest Ergoweale Flas Left 30 cm           | A11202250   |
| 9   | Seat ErgoMedic Plus 29x32 cm           | ST205546    |         |   |             |
| 9   | Seat ErgoMedic Plus 33x35 cm           | ST205046    |         |   |             |
| 9   | Seat ErgoMedic Plus 36x39 cm           | ST203046    |         |   |             |
| 9   | Seat ErgoMedic Plus 36x43 cm           | ST204046    |         |   |             |
| 9   | Seat ErgoMedic Plus 30x43 cm           | ST204046    |         | .1  |             |
| 9   | Jear Ligorieule i lus 40x33 cm         | 01200040    |         | ★ The backrests are also avaliable in Slowr | ecovery.    |



#### OTHER UPHOLSTERED PARTS

#### Trunk support

The trunk support is used to produce more upright, active sitting. Mercado has three types: Curved and adjustable, straight, and curved trunk support.

#### Neck rest

The neckrest provides support for the head in order to achieve the best sitting position and to provide relief for the neck and spine even for users with poor body control. It can easily be adjusted to achieve the best setting for the individual user. Available in flat, curved or ERGOMEDIC® PLUS versions.

#### Cover

The covers protect the upholstery from dirt and dust, as well as extending the upholstery's lifetime. Covers are available in different colours for a variety of combinations.

Incontinency covers are a quick and easy solution to protect seats against moisture.

#### Armrest

The Hemiplegi armrest is suitable for partially paralysed, or users with restricted control. The armrest panel has been hollowed out, which makes it easier to keep the paralysed arm on the rest. The armrest panel also rotates, which allows for the best possible position. The mobile armrest panel is upholstered with extra foam for optimised comfort.

#### Side support

The side supports help the user maintain stability and position. These can also be used as thigh supports, so that the users feet are better supported on the footrests or the floor, and prevents their legs from sliding apart. It can be adjusted for height and width, or turned backwards and forwards. Available in three sizes: S, M and L.

#### Slow Recovery seat and others

Slow Recovery seats are manufactured with two materials in the cushioning, which provides optimum support and comfort to the user. The upper layer comprises of a memory foam that moulds itself, and provides good pressure distribution and increased sitting comfort, while the lower layer raises the seat bone nodules from the wood plate. Even if the user does not sit in the chair for a sufficiently long period of time, the seat cushion adjusts itself so that a cavity is formed where the person sits. This material is available for all upholstered parts.

The wooden seat disc should be used when casting or retrofitting of its shell. The seat with anti-slip is recommended when using individually adjusted seats. Seat with hygiene fabric should be used by users with incontinence problems. To be completed with a cover in plush.



| Trunk support, width x height       |                                     |               | Article No.     | Ne  | ckrest, width x h   | _                 | Article No. |
|-------------------------------------|-------------------------------------|---------------|-----------------|---|---------------------|-------------------|-------------|
| -                                   | 1 Trunk support cushion 1           |               | TS1046          | 3   |                     | Straight 29x16 cm |             |
| 2                                   | 2 Trunk support curved and          | d adjustable  | TS2038          | 4   | Neckrest cushion    | Curved 31x18 cm   | NS3046      |
| 1                                   | Armrest                             |               | Article No.     | Art   | icle No. for angle  | •                 |             |
| 6                                   | 6 Armrest ErgoMedic 25 c            | m             | AR101038        | 240387  |                     |                   |             |
| 6                                   | 6 Armrest ErgoMedic 30 c            | m             | AR102038        | 240386  |                     |                   |             |
| 6                                   | 6 Armrest ErgoMedic 35 c            | m             | AR103038        | 240385  |                     |                   |             |
| 6                                   | 6 Armrest ErgoMedic 40 c            | m             | AR104038        | 240388  |                     |                   |             |
| 7                                   | 7 Armrest Mobile 25 cm              |               | AR301038        | 245387  |                     |                   |             |
| 7                                   | 7 Armrest Mobile 30 cm              |               | AR302038        | 245386  |                     |                   |             |
| 7                                   | 7 Armrest Mobile 35 cm              |               | AR303038        | 24  | 245385              |                   |             |
| 7                                   | 7 Armrest Mobile 40 cm              |               | AR304038        | 245388  |                     |                   |             |
| 8 Armrest Hemiplegic 45x12 cm       |                                     | AR406038      | 244387          |   |                     |                   |             |
| ę                                   | 9 Wooden seat disc                  |               | ST10xx90        | Covers  |                     |                   |             |
| 10 Wooden seat disc with anti-slip  |                                     | ST10xx61      | 5               | Exist in plush, patterned fabric and incontinence |                     |                   |             |
| 9                                   | Side support, width x height        |               | Article No.     | material. Log on to www.mercado-medic             |                     | medic.co.uk or    |             |
| 11 Side support cushion XS, 15x8 cm |                                     | SS1046        |                 | material. contac                                  | ct Mercado for more | e information.    |             |
| -                                   | 12 Side support cushion S, 28x8 cm  |               | SS2046          |   |                     |                   |             |
| -                                   | 13 Side support cushion M, 24x14 cm |               | SS3046          | 16 Seat ErgoMedic with hygiene                    |                     |                   |             |
| -                                   | 14 Side support cushion L,          | 30x16 cm      | SS4046          |   | completed with      | cover in plush    |             |
| 9                                   | Seat Slow recovery                  | Slow recovery | Slow recovery   | Slo   | w recovery          | Slow recovery     |             |
| 9                                   | Size in width x depth               | ErgoMedic     | ErgoMedic coxit | Erg   | oMedic PLUS         | ErgoMedic PLUS    | coxit       |
| -                                   | 15 33x35 cm                         | ST115046      | -               | ST  | 215046              | -                 |             |
| -                                   | 15 36x39 cm                         | ST114546      | -               | ST  | 214546              | -                 |             |
| •                                   | 15 36x43 cm                         | ST114046      | -               | ST  | 214046              | -                 |             |
| •                                   | 15 40x39 cm                         | ST113546      | ST143546        | ST  | 213546              | ST243546          |             |
| -                                   | 15 40x43 cm                         | ST113046      | ST143046        | ST  | 213046              | ST243046          |             |
|                                     |                                     |               |                 |   |                     |                   |             |

ST142546

ST142046

ST141546

ST141046

ST212546

ST212046

ST211546

ST211046

ST242546

ST242046

ST241546

ST241046

15

40x48 cm

15 44x43 cm

15 44x48 cm

15 48x48 cm

ST112546


ST112046


ST111546

ST111046

#### MEASUREMENTS FOR BACKRESTS

#### ErgoMedic





 1
 Size: 42x43 cm
 Art.nr: BR101546

 2
 Size: 40x34 cm
 Art.nr: BR103046


 3
 Size: 38x26 cm
 Art.nr: BR105546


 4
 Size: 33x36 cm
 Art.nr: BR102546

 5
 Size: 33x26 cm
 Art.nr: BR106046

6 Size: 33x16 cm Art.nr: BR109046 7 Size: 30x34 cm Art.nr: BR103546 8 Size: 30x26 cm Art.nr: BR106546 9 Size: 27x32 cm Art.nr: BR104546 10 Size: 27x22 cm Art.nr: BR107046


#### ErgoMedic PLUS






1 Size: 46x43 cm Art.nr: BR201046 2 Size: 42x43 cm Art.nr: BR201546 3 Size: 35x43 cm Art.nr: BR202046 4 Size: 40x20 cm Art.nr: BR207546 5 Size: 35x20 cm Art.nr: BR208046 6 Size: 30x34 cm Art.nr: BR203546 7 Size: 25x34 cm Art.nr: BR204046 8 Size: 23x27 cm Art.nr: BR205046 9 Size: 23x20 cm Art.nr: BR208546

#### ErgoMedic PLUS with wedge



#### ErgoMedic PLUS w reinforced wedge



 1 Size: 46x43 cm
 Art.nr: BR251046

 2 Size: 42x43 cm
 Art.nr: BR251546

 3 Size: 35x43 cm
 Art.nr: BR252046


 4 Size: 30x34 cm
 Art.nr: BR253546

1 Size: 46x43 cm Art.nr: BR351046 2 Size: 42x43 cm Art.nr: BR351546 3 Size: 35x43 cm Art.nr: BR352046

#### ErgoMedic PLUS Comfort

# 43 cm 24 cm 43 cm 43 cm 35 cm


#### ErgoMedic PLUS Comfort with wedge




1 Size: 46x43 cm Art.nr: BR211046 2 Size: 42x43 cm Art.nr: BR211546 3 Size: 35x43 cm Art.nr: BR212046 1 Size: 46x43 cm Art.nr: BR261046 2 Size: 42x43 cm Art.nr: BR261546 3 Size: 35x43 cm Art.nr: BR262046

#### **MEASUREMENTS FOR SEATS**


#### ErgoMedic






1 Size: 48x48 cm Art.nr: ST101046 2 Size: 44x48 cm Art.nr: ST101546 3 Size: 44x43 cm Art.nr: ST102046 4 Size: 40x48 cm Art.nr: ST102546 5 Size: 40x43 cm Art.nr: ST103046 6 Size: 40x39 cm Art.nr: ST103546 7 Size: 36x43 cm Art.nr: ST104046 8 Size: 36x39 cm Art.nr: ST104546 9 Size: 33x35 cm Art.nr: ST105046 10 Size: 29x32 cm Art.nr: ST105546

#### ErgoMedic Slowrecovery






1 Size: 48x48 cm Art.nr: ST111046 2 Size: 44x48 cm Art.nr: ST111546 3 Size: 44x43 cm Art.nr: ST112046 4 Size: 40x48 cm Art.nr: ST112546 5 Size: 40x43 cm Art.nr: ST113046

6 Size: 40x39 cm Art.nr: ST113546 7 Size: 36x43 cm Art.nr: ST114046 8 Size: 36x39 cm Art.nr: ST114546 9 Size: 33x35 cm Art.nr: ST115046

#### ErgoMedic Coxit

# 1 48 cm 48 cm 48 cm 43 cm 43 cm

#### ErgoMedic Coxit Slowrecovery



| 1 | Size: 48x48 cm | Art.nr: ST131046 |
|---|----------------|------------------|
| 2 | Size: 44x48 cm | Art.nr: ST131546 |
| 3 | Size: 44x43 cm | Art.nr: ST132046 |
| 4 | Size: 40x48 cm | Art.nr: ST132546 |
| 5 | Size: 40x43 cm | Art.nr: ST133046 |
| 6 | Size: 40x39 cm | Art.nr: ST133546 |


40 cm ~


39 cm

43 cm

39 cm

#### ErgoMedic PLUS







1 Size: 48x48 cm Art.nr: ST201046 2 Size: 44x48 cm Art.nr: ST201546 3 Size: 44x43 cm Art.nr: ST202046 4 Size: 40x48 cm Art.nr: ST202546 5 Size: 40x43 cm Art.nr: ST203046

6 Size: 40x39 cm Art.nr: ST203546 7 Size: 36x43 cm Art.nr: ST204046 8 Size: 36x39 cm Art.nr: ST204546 9 Size: 33x35 cm Art.nr: ST205046 10 Size: 29x32 cm Art.nr: ST205546

#### ErgoMedic PLUS Slowrecovery






1 Size: 48x48 cm Art.nr: ST211046 2 Size: 44x48 cm Art.nr: ST211546 3 Size: 44x43 cm Art.nr: ST212046 4 Size: 40x48 cm Art.nr: ST212546 5 Size: 40x43 cm Art.nr: ST213046 6 Size: 40x39 cm Art.nr: ST213546 7 Size: 36x43 cm Art.nr: ST214046 8 Size: 36x39 cm Art.nr: ST214546 9 Size: 33x35 cm Art.nr: ST215046

#### ErgoMedic PLUS Coxit


## 21 cm 48 cm 12 cm 48 cm 17 cm 48 cm 48 cm 13 cm 48 cm 43 cm 43 cm 13 cm 43 cm 43 cm 13 cm 39 cm

#### ErgoMedic PLUS Coxit Slowrecovery




1 Size: 48x48 cm Art.nr: ST231046 2 Size: 44x48 cm Art.nr: ST231546 3 Size: 44x43 cm Art.nr: ST232046 4 Size: 40x48 cm Art.nr: ST232546 5 Size: 40x43 cm Art.nr: ST233046 6 Size: 40x39 cm Art.nr: ST233546 1 Size: 48x48 cm Art.nr: ST241046 2 Size: 44x48 cm Art.nr: ST241546 3 Size: 44x43 cm Art.nr: ST242046 4 Size: 40x48 cm Art.nr: ST242546 5 Size: 40x43 cm Art.nr: ST243046 6 Size: 40x39 cm Art.nr: ST243546

#### Wooden seat disc





#### Seat with anti-slip




1 Size: 48x48 cm Art.nr: ST101090 2 Size: 44x48 cm Art.nr: ST101590 Size: 44x43 cm Art.nr: ST102090 4 Size: 40x48 cm Art.nr: ST102590 5 Size: 40x43 cm Art.nr: ST103090 Size: 40x39 cm Art.nr: ST103590 Size: 36x43 cm Art.nr: ST104090 8 Size: 36x39 cm Art.nr: ST104590 Size: 33x35 cm Art.nr: ST105090

1 Size: 48x48 cm Art.nr: ST101061 2 Size: 44x48 cm Art.nr: ST101561 Size: 44x43 cm Art.nr: ST102061 4 Size: 40x48 cm Art.nr: ST102561 5 Size: 40x43 cm Art.nr: ST103061 Size: 40x39 cm Art.nr: ST103561 Size: 36x43 cm Art.nr: ST104061 8 Size: 36x39 cm Art.nr: ST104561 Size: 33x35 cm Art.nr: ST105061











#### EXPLANATION OF THE STRUCTURE OF THE ARTICLE NUMBER Art.no: ST101546 = "Seat ErgoMedic, plush grey 44x48 cm"

Digit 1-2 = Type of upholstery article(ST = seat, BR = back rest etc.).

Digit 3-4 = Seatsystem (ErgoMedic, ErgoMedic PLUS, etc.).

Digit 5-6 = Size

Digit 7 = Material (plush, skai etc).

Digit 8 = Colour

#### **UPHOLSTERY**

Plush

41 Bordeaux 42 Brown 44 Marine blue 45 REAL-blue 46 Grey 48 Black

Skai

31 Red 32 Brown 35 Blue 36 Grey 38 Black

Patterned fabric, Adult Patterned fabric, Child

75 Blue 81 Red 85 Blue 89 Yellow

Alternative covering

60 Hygiene fabric 95 Without clothing

Covering

66 Incontinence cover, grey plastic-coated towel 73 Child/Yellow cotton 74 Child/Blue cotton

Coverings are available for all of Mercado Medic's upholstered parts, visit www.mercado-medic.co.uk or contact us for more information.

