

**TOPRO
TAURUS**

TOPRO



STEADY!

Safe walkers for institutions and at home

www.topro.co.uk

Steady and safe walkers for institutions and at home

- Easy to read heightscale.
- Easy to adjust the height.
- High enough for tall users (up to 2 m).
- Footbrake on all four wheels.
- Tested and approved according to ISO 11199-3 for use indoors on plain foundation.
- The TOPRO TAURUS E-models are tested and approved according to IEC 60601-1:2005 + corr 1:2006 + corr 2:2007 and IEC 60601-1-2:2007

Maintenance free material

- Aluminium with polyester surface and polyurethane.
- Easy to clean.
- Non-corrosive.
- All wear and tear parts can be exchanged.

Two models

TOPRO TAURUS H Hydraulic height adjustment.

TOPRO TAURUS E Electrical height adjustment.

The two models are available in three different types:



H-JHB / E-JHB	YES	YES	YES	YES
H-JH / E-JH	YES	-	-	-
H / E	-	-	-	-

The models without brakes are delivered partly assembled.

	TOPRO ART. No.
TOPRO TAURUS H-JHB	814787
TOPRO TAURUS E-JHB	814789
TOPRO TAURUS H-JH	814786
TOPRO TAURUS E-JH	814788
TOPRO TAURUS H	814785
TOPRO TAURUS E	814790

TOPRO TAURUS H

- Gas springs simplify the height adjustment.
- Can be used in the shower.

Product weight 18 kg / 2 ¾ Stone

Max user weight 150 kg / 23 ½ Stone





Max 22 cm / 8 3/4\"/>

66 cm / 26\"/>

78 cm / 30 3/4\"/>

91 cm / 35 3/4\"/>



34 cm / 13 1/2\"/>

34 - 45 cm / 13 1/2\"/>

3,6 cm / 1 1/2\"/>



12.5 cm / 5 1/8\"/>

3,1 cm / 1 1/4\"/>

**TOPRO
TAURUS E**

- Electric motor makes the height adjustment very simple.
- Can also be used as raising support.

Product weight 22 kg / 3 1/2 Stone

Max user weight 150 kg / 23 1/2 Stone

Accessories for a variety of uses

To be ordered separately.

Product **TOPRO ART. No.**



Tray **814810**



Foot plate **814645**



Oxygen bottle holder **814820**



Brake handle left side **814846**



Stand for body fluid-boxes **814034**



Height adjustment panel **814805**



Adjustable handles **814001**



Small basket **814821**



Infusion stand **814793**